LABORATORY OF FOREST PRODUCT CHEMISTRY FACULTY OF FORESTRY MULAWARMAN UNIVERSITY

PROTOCOL 009

STANDARD OPERATING PROCEDURE

	Date:
Title: <i>Total Phenolic Content</i> Testing Method	Created by:
	File:

1. OBJECTIVE

To provide instructions for performing a total phenolic content test.

2. PRINCIPLE

The total phenolic content test is conducted to determine the total phenolic content in a sample using a colorimetric method.

3. REQUIRED EQUIPMENT

Laboratory suits, safety glasses, clean wipes/tissues.

3. MATERIALS AND EQUIPMENT

Material	EQUIPMENT	
Plant extract sample	Test tube	
Phenol standards (Example: Galic acid) Analytical balances		
Folin-Ciocalteu Solution	UV-Vis Spectrophotometer	
Sodium carbonate (Na2CO3)	Beaker glass	
Distilled water	Measuring cup	
	Micropipettes & Tips	

4. PROCEDURE

4.1. PREPARATION OF TEST SOLUTION

- 4.1.1 Weigh 1 mg of extract, dissolve in 10 ml of distilled water
- 4.1.2 Weigh 1 mg of standard phenol, dissolve in 10 ml of distilled water
- 4.1.3 Dissolve 1 ml of Folin-Ciocalteu solution in 10 ml of distilled water
- 4.1.4 Weigh 7.5 g of Na2CO3, dissolve in 100 ml of distilled water

4.2. TESTING OF TOTAL PHENOLIC CONTENT

4.2.1 To prepare the phenolic standard curve, add the solutions to each test tube in the order specified in the table below;

Tube	Phenol Standard	Distilled water (µl)	Folin-Ciocalteu (µl)	Na2CO3 (µl)
	(μl)			
Tube 1	0	500	250	1250
Tube 2	20	480	250	1250
Tube 3	40	460	250	1250
Tube4	60	440	250	1250
Tube 5	80	420	250	1250
Tube 6	100	400	250	1250

4.2.2 For the test sample, add the solutions to each test tube as per the table below, and repeat this step three times;

Tube	Test Sample (µl)	Distilled water (µl)	Folin-Ciocalteu (µl)	Na2CO3 (µl)
Tube 1	100	400	250	1250

4.2.3 Incubate the solutions for 60 minutes.

4.2.4 Absorbance with a wavelength of 760 nm

4.2.5 Create a calibration curve from the absorbance results of the phenolic standards to obtain values and a regression equation.;

4.2.6 Insert the average absorbance value of the sample into the regression equation;

4.2.7 Calculate the total phenol content using the formula;

$$C = c \times \frac{v}{w}$$

Where: C : Total phenol content (μg (standard phenol)/mg extract)

- c : Concentration of sample taken (µg/ml)
- m : Extract sample weight (μg)
- v : Volume (ml)

4.2.8 Clean the work area and equipment if you have finished testing.

4.3. CAUTION

- 4.3.1 Weigh all solutions and conduct testing in a low-light area to prevent rapid reactions in the solutions;
- 4.3.2 It is advisable to weigh the solutions one day before testing to stabilize reactions occurring due to dilution;4.3.3 Store solutions in a refrigerator;

Version : 06 - 2023 Validated by :